



# MYT-3CRP CAL

HEMATOLOGY CALIBRATOR

**CALIBRATOR**

**LOT** N0524

 2024-06-10

## ASSAY VALUES

Instrument	Parameter	Unit	Reference Value
Orphée Mythic 3CRP	WBC	$\times 10^3/\mu\text{L}$	9.78
	RBC	$\times 10^6/\mu\text{L}$	4.59
	HGB	g/dL	13.3
	MCV	fL	89.1
	PLT	$\times 10^3/\mu\text{L}$	243

**NOTE** Before using, refer to the instruction sheet for mixing directions.  
All brands and products are trademarks or registered trademarks of their respective companies.

**ORPHÉE SA**  
19 chemin du Champ-des-Filles  
CH-1228 Plan-les-Ouates / Geneva  
Switzerland

**EC REP**

**PZ CORMAY S.A.**  
ul. Wiosenna 22,  
05-092 Łomianki, POLSKA

**CE** **IVD**