



MYT-3CRP

HEMATOLOGY CONTROLS

CONTROL

LOT B0524



2024-08-10

ASSAY VALUES AND EXPECTED RANGES

Instrument	Parameter	Low		Normal		High	
		LOT	B0524L	LOT	B0524N	LOT	B0524H
Orphée Mythic 3CRP	WBC × 10 ³ /μL		2.14 ± 0.60		8.01 ± 1.00		21.84 ± 2.50
	LYM %		64.3 ± 9.0		36.0 ± 8.0		16.7 ± 8.0
	GRAN %		32.1 ± 9.0		60.4 ± 8.0		72.5 ± 8.0
	MID %		3.6 ± 3.6		3.6 ± 3.6		10.8 ± 6.0
	LYM × 10 ³ /μL		1.37 ± 0.30		2.88 ± 0.70		3.65 ± 2.00
	GRAN × 10 ³ /μL		0.69 ± 0.30		4.84 ± 0.70		15.83 ± 2.00
	MID × 10 ³ /μL		0.08 ± 0.08		0.29 ± 0.29		2.36 ± 1.00
	RBC × 10 ⁶ /μL		2.39 ± 0.24		4.65 ± 0.24		5.79 ± 0.30
	HGB g/dL		5.9 ± 0.6		13.5 ± 0.6		18.9 ± 0.8
	HCT %		19.0 ± 2.0		41.4 ± 3.0		55.9 ± 4.0
	MCV fL		79.6 ± 5.0		89.1 ± 5.0		96.5 ± 5.0
	MCH pg		24.1 ± 2.5		28.9 ± 2.5		32.4 ± 2.5
	MCHC g/dL		31.3 ± 3.0		33.2 ± 3.0		34.6 ± 3.0
	RDW-CV %		17.0 ± 3.0		15.1 ± 3.0		14.3 ± 3.0
	RDW-SD fL		49.6 ± 8.0		49.4 ± 8.0		52.4 ± 8.0
	PLT × 10 ³ /μL		74 ± 20		268 ± 40		516 ± 60
	MPV fL		9.3 ± 3.0		9.5 ± 3.0		9.4 ± 3.0
	PDW fL		12.0 ± 3.0		12.6 ± 3.0		12.6 ± 3.0
	PCT %		0.069 ± 0.069		0.256 ± 0.200		0.491 ± 0.200
	P-LCR %		26.5 ± 10.0		28.3 ± 10.0		28.1 ± 10.0
P-LCC × 10 ³ /μL		20 ± 15		77 ± 25		146 ± 30	

NOTE Before using, refer to the instruction sheet for mixing directions.
All brands and products are trademarks or registered trademarks of their respective companies.



B0524L



B0524N



B0524H

ORPHÉE SA
19 chemin du Champ-des-Filles
CH-1228 Plan-les-Ouates / Geneva
Switzerland



PZ CORMAY S.A.
ul. Wiosenna 22,
05-092 Łomianki,
POLSKA

