




# MYTHIC CRP Controls

CRP CONTROLS

**CONTROL**

**LOT** 2024032101

 2025-03-20

## ASSAY VALUES AND EXPECTED RANGES

Instrument	Control	Reference Value (mg/L)	Deviation Range (mg/L)
Orphée Mythic 3CRP	Level 1	3.4	± 1.1
	Level 2	14.6	± 5.1
	Level 3	41.3	± 12.7

**NOTE** Before using, refer to the instruction sheet for mixing directions.  
All brands and products are trademarks or registered trademarks of their respective companies.

**ORPHÉE SA**  
19 chemin du Champ-des-Filles  
CH-1228 Plan-les-Ouates / Geneva  
Switzerland

**EC REP**

**PZ CORMAY S.A.**  
ul. Wiosenna 22,  
05-092 Łomianki, POLSKA

**CE** **IVD**