



Myt-5D VET

HEMATOLOGY CONTROLS

CONTROL

MYT503-VET
LOT MYTV2405
 2024/07/05

ASSAY VALUES AND EXPECTED RANGES

Instrument	Parameter	Low	Normal	High
		LOT MYTV2405L	LOT MYTV2405N	LOT MYTV2405H
Orphée Mythic 5VET PRO	WBC (WOC) × 10 ⁹ /L	2.73 ± 0.4	6.86 ± 0.8	17.05 ± 2.2
	NEU #	1.35 ± 0.60	4.46 ± 1.20	12.45 ± 3.40
	NEU %	50.95 ± 10.00	65.07 ± 8.00	72.99 ± 8.00
	LYM #	1.03 ± 0.50	1.56 ± 0.90	2.26 ± 1.70
	LYM %	37.79 ± 9.00	22.77 ± 7.00	13.25 ± 6.00
	MON #	0.22 ± 0.22	0.50 ± 0.50	0.91 ± 0.91
	MON %	7.96 ± 6.40	7.26 ± 6.50	5.32 ± 5.32
	EOS #	0.11 ± 0.11	0.31 ± 0.31	1.38 ± 1.38
	EOS %	2.80 ± 2.80	4.58 ± 4.58	8.08 ± 7.30
	BASO #	0.02 ± 0.02	0.02 ± 0.02	0.06 ± 0.06
	BASO %	0.51 ± 0.51	0.32 ± 0.32	0.35 ± 0.35
	RBC × 10 ¹² /L	1.90 ± 0.15	4.35 ± 0.20	4.95 ± 0.24
	HGB g/dL	5.7 ± 0.5	14.2 ± 0.5	16.8 ± 1.0
	HGB g/L	57.0 ± 5	142.0 ± 5	167.8 ± 10
	HCT %	14.9 ± 1.8	39.0 ± 2.4	45.8 ± 3.0
	MCV fL	78.6 ± 4.0	89.7 ± 4.0	93.6 ± 4.0
	MCH pg	30.2 ± 2.8	32.7 ± 2.0	34.3 ± 2.0
	MCHC g/dL	38.3 ± 3.6	36.5 ± 2.8	36.6 ± 2.8
	MCHC g/L	383 ± 36	365 ± 28	366 ± 28
	RDW-CV %	17.8 ± 3.0	15.1 ± 3.0	14.7 ± 3.0
PLT × 10 ⁹ /L	80 ± 20	229 ± 40	425 ± 65	
MPV fL	8.7 ± 3.0	8.7 ± 3.0	8.8 ± 3.0	

ORPHEE SA
 19 chemin du Champ-des-Filles
 CH-1228 Plan-les-Ouates / Geneva
 Switzerland