

Lublin, Poland, 13th February 2024**WITHDRAWAL FROM THE OFFER**

Dear Valuable Customers,

We would like to inform that due to offer optimization, it has been made decision about withdrawal of **triglycerides mono** reagent (**TG mono**).

Until now, PZ CORMAY S.A. had in the offer two types of reagents for determination of triglycerides concentration used both for manual assay and in several automatic analysers, as below:

- ✓ TG mono-reagent,
- ✓ TG double-reagent.

Principle of the method for both parameters is based on similar colorimetric, enzymatic method with glycerophosphate oxidase.

The decision to leave in offer just one, more popular TG double-reagent has been made. TG mono will be available until stock last. Then, we recommend to order TG double-reagent instead. Please, check Table 1 and 2 with more details about discontinued and potential replacement kits.

Tab 1. TG mono to be withdrawn from the offer

Cat. No	Description
2-282	Liquick Cor-TG mono 30
2-283	Liquick Cor-TG mono 60
2-284	Liquick Cor-TG mono 120
7-273	ACCENT-200 TG mono
7-373	ACCENT-300 TG mono
7-473	A-400 TG mono
4-226	PRESTIGE 24i LQ TG mono
4-426	PRESTIGE 24i LQ TG mono
4-573	HC-TG mono
9-473	OS-TG mono
5-518	B50-TG mono
7-810	A-800 TG mono
2-286_R1	Liquick Cor-TG mono BULK

Tab 2. Potential replacement

Cat. No	Description
2-262	Liquick Cor-TG 30
2-253	Liquick Cor-TG 60
2-254	Liquick Cor-TG 120
7-253	ACCENT-200 TG
7-353	ACCENT-300 TG
7-453	A-400 TG
4-253	PRESTIGE 24i LQ TG
4-453	PRESTIGE 24i LQ TG
4-553	HC-TG
9-407	OS-TG
5-517	B50-TG
7-809	A-800 TG
2-280_R1	Liquick Cor-TG BULK
2-280_R2	Liquick Cor-TG BULK

We apologize for any inconvenience.

For any further information, please contact your dedicated Area Sales Manager.

Yours sincerely,

Clinical Chemistry Product Manager

W/1/2024

PZ Cormay S.A.

Poland, 05-092 Łomianki 22 Wiosenna Street
NIP: 1181872269
REGON: 140777556

Correspondence address:

PZ Cormay S.A.
Poland, 02-785 Warsaw, 303 Pulawska Street

office@cormay.com
tel. +48 (22) 751 79 10